

GenJet™ In Vitro DNA Transfection Reagent (Ver. II)

----- A General Protocol for Transfecting Mammalian Cells

- 100 µl
 500 µl
 1000 µl



10075 Tyler Place, Suite 19
 Ijamsville, MD 21754
 FAX. 301-560-4919
 TEL. 301-330-5966
 Toll Free. 1-(866)-918-6812
 Email: info@signagen.com
 Web: www.signagen.com

This product is for laboratory research ONLY and not for diagnostic use

Introduction:

GenJet™ In Vitro DNA Transfection Reagent (Ver. II) is upgraded version of GenJet™ In Vitro DNA Transfection Reagent. With a new chemistry, more DNA condensing groups were released in the new version compared with old version GenJet™, leading to 3-4 times more efficient in DNA delivery. GenJet™ (Ver. II) was shown to deliver genes to various established cell lines as well as primary cells.

Important Transfection Guidelines for New Version:

- Do NOT follow transfection procedures for GenJet old version. Read protocol for new version carefully before transfection
- For high efficiency, transfect cells at high density. 70-80% confluency is highly recommended
- To lower cytotoxicity, transfect cells in presence of serum (10%) and antibiotics
- An advanced protocol is provided for hard-to-transfect cells only if general protocol gives less than 10% efficiency (see back page)

Part I. General Protocol for Transfecting Mammalian Cells:

Step I. Cell Seeding:

Cells should be plated 18 to 24 hours prior to transfection so that the monolayer cell density reaches to the optimal 70-80% confluency at the time of transfection. Complete culture medium with serum and antibiotics is freshly added to each well 30-60 minutes before transfection.

Note: High serum levels (>5%) with antibiotics usually do not have inhibitory effect on transfection efficiency. For some specific cells, maximal transfection efficiencies are observed in the presence of serum and antibiotics. We recommend using complete serum/antibiotics-containing medium initially.

Table 1. Recommended Amounts for Different Culture Vessel Formats

| Culture Dish | Culture Medium (ml) | Plasmid DNA (µg) | Diluent Volume (mL) | GenJet™ Reagent (µL) |
|---------------|---------------------|------------------|---------------------|----------------------|
| 48 well plate | 0.3 | 0.25 | 2 x 0.015 | 0.75 |
| 12 well plate | 0.75 | 0.75 | 2 x 0.038 | 2.25 |
| 6-well plate | 1.0 | 1 | 2 x 0.05 | 3.0 |
| 35 mm dish | 1.0 | 1 | 2 x 0.05 | 3.0 |
| 60 mm dish | 2.8 | 2.5 | 2 x 0.10 | 7.5 |
| 10 cm dish | 5.0 | 5 | 2 x 0.25 | 15 |
| T75 flask | 8.0 | 9 - 18 | 2 x 0.40 | 27 - 54 |
| 250 ml flask | 18 | 25 - 50 | 2 x 0.8 | 75 - 150 |

Step II. Preparation of GenJet™-DNA Complex and Transfection Procedures

For different cell types, the optimal ratio of GenJet™ (µL):DNA (µg) is around 3:1. We recommend the GenJet™ (µL):DNA (µg) ratio of 3:1 as a starting point which usually gives satisfactory transfection efficiency with

invisible cytotoxicity. To ensure the optimal size of GenJet™/DNA complex particles, we recommend using serum-free DMEM with High Glucose to dilute DNA and GenJet™ Reagent.

The following protocol is given for transfection in 24-well plates, refer to **Table 1** for transfection in other culture formats. The optimal transfection conditions for a majority of adherent cell lines are given in the standard protocol described below.

- For each well, add 0.5 ml of complete medium with serum and antibiotics freshly 30-60 minutes before transfection.
- For each well, dilute 0.5 µg of DNA into 25 µl of serum-free DMEM with High Glucose. Vortex gently and spin down briefly to bring drops to bottom of the tube .
- For each well, dilute 1.5 µl of GenJet™ reagent (Ver. II) into 25 µl of serum-free DMEM with High Glucose. Pipette up and down 3-4 times to mix.
- Note:** Never use Opti-MEM to dilute GenJet™ reagent and DNA, it will disrupt transfection complex.
- Add the diluted GenJet™ Reagent immediately to the diluted DNA solution all at once. (**Important: do not mix the solutions in the reverse order !**)
- Immediately pipette up and down 3-4 times or vortex to mix followed by incubation for-10 minutes at room temperature to allow GenJet™-DNA complexes to form.
- Note:** Never keep the DNA/GenJet™ complex longer than 20 minutes
- Add the 50 µl GenJet™/ DNA complex drop-wise onto the medium in each well and homogenize the mixture by gently swirling the plate.
- Remove DNA/GenJet™ complex and replace with complete serum/antibiotics containing medium 12-18 hours post transfection. **For sensitive cells, to lower cytotoxicity, remove GenJet™/DNA complex and replace with complete medium 5 hours after transfection.**
- Check transfection efficiency 24 to 48 hours post transfection.

Storage: GenJet™ DNA In Vitro Transfection Reagent is stable for up to 12 months at +4 °C. This item shipped at ambient temperature

GenJet™ DNA In Vitro Transfection Reagent (Ver. II)

----- An Advanced Protocol for Transfecting Hard-to-Transfect Mammalian Cells

- 100 µl
- 500 µl
- 1000 µl



10075 Tyler Place, Suite 19
Ijamsville, MD 21754
FAX. 301-560-4919
TEL. 301-330-5966
Toll Free. 1-(866)-918-6812
Email: info@signagen.com
Web: www.signagen.com

This product is for laboratory research ONLY and not for diagnostic use

Part II. Advanced Protocol for Transfecting Hard-To-Transfect Mammalian Cells

Important: The advanced protocol for hard-to-transfect cells is provided only if general protocol gives less than 10% efficiency. **For some primary cells which cannot be trypsinized (like primary neurons), go directly to Step II, skip trypsinization and incubate freshly prepared primary cell pellet with transfection complex.**

Step I. Cell Culture Before Transfection

Cells should be plated at least 24 hours prior to transfection so that the monolayer cell density reaches to the optimal 95-100% confluency at the day of transfection.

Table 2. A Guideline for Optimal Cell Number Per Well in Different Culture Formats

| Culture Dishes | Surface Area (cm ²) | Optimal Cell Number |
|----------------|---------------------------------|------------------------|
| T75 Flask | 75 | 9.6 x 10 ⁶ |
| 100 mm Dish | 58 | 7.3 x 10 ⁶ |
| 60 mm Dish | 21 | 2.7 x 10 ⁶ |
| 35 mm Dish | 9.6 | 1.2 x 10 ⁶ |
| 6-well Plate | 9.6 | 1.2 x 10 ⁶ |
| 12-well Plate | 3.5 | 0.44 x 10 ⁶ |
| 24-well Plate | 1.9 | 0.24 x 10 ⁶ |
| 48-well Plate | 1.0 | 0.11 x 10 ⁶ |
| 96-well Plate | 0.3 | 0.31 x 10 ⁵ |

Table 3. Recommended Amounts for Different Culture Vessel Formats

| Culture Dish | Transfection Complex Volume (ml) | Plasmid DNA (µg) | GenJet™ Reagent (µL) |
|--------------|----------------------------------|------------------|----------------------|
| 96-well | 0.02 | 0.2 | 0.8 |
| 48-well | 0.04 | 0.5 | 2 |
| 24-well | 0.1 | 1.0 | 4 |
| 6-well | 0.2 | 2 | 8 |
| 35 mm dish | 0.2 | 2 | 8 |
| 60 mm dish | 0.5 | 5 | 20 |
| 10 cm dish | 1.0 | 8 | 32 |
| T75 flask | 1.5 | 36 | 144 |

Step II. Preparation of Cells in Suspension

The following protocol is given for transfecting hard-to-transfect mammalian cells in 6-well plates, refer to [Table 2](#) for optimal cell

number per well per culture vessels' surface area. The optimal transfection conditions for mammalian cells are given in the standard protocol described below.

- Detach the cells with trypsin/EDTA and stop the trypsinization with complete culture medium.

Note: Cells that are difficult to detach may be placed at 37 °C for 5-15 min to facilitate detachment

- Take an aliquot of trypsinized cell suspension and count the cells to determine the cell density.

- Centrifuge the required 1.2x10⁶ cells per well for 6-well plate at 150xg at room temperature for 10 min.

- Use fine tip pipette to remove supernatant **completely** so that no residual medium covers the cell pellet.

Step III. Preparation and Application of Transfection Complex

For hard-to-transfect mammalian cells, the optimal ratio of GenJet™ (µL):DNA (µg) is 4:1. To ensure the optimal size of complex particles, we recommend using serum-free DMEM with High Glucose to dilute DNA and GenJet™ Reagent.

The following protocol is given for transfection in 6-well plates, refer to [Table 3](#) for transfection in other culture formats.

- For each well of 6-well plate, dilute 2 µg of DNA into 100 µl of serum-free DMEM with High Glucose. Vortex gently and spin down briefly to bring drops to bottom of the tube.

- For each well of 6-well plate, dilute 8 µl of GenJet™ reagent (Ver. II) into 100 µl of serum-free DMEM with High Glucose. Vortex gently and spin down briefly.

- Add the diluted GenJet™ Reagent immediately to the diluted DNA solution all at once.

- Immediately pipette up and down 3-4 times or vortex briefly to mix followed by incubation for ~15 minutes at room temperature to allow transfection complexes to form.

Important: Never keep the transfection complexes longer than 15 minutes

- **Gently** resuspend the cell pellet prepared from [Step II](#) immediately in the 200 µl transfection complex and incubate at 37 °C for 20 minutes.

- At the end of incubation, add 2.0 ml of pre-warmed fresh complete cell growth medium to cells and plate onto one well of a 6-well plate.

- Remove transfection complex containing medium gently and refill with complete culture medium 5-12 hours after plating.

- Check efficiency 24 to 48 hours post transfection.