

PowerFect™ In Vitro DNA and siRNA Transfection Kit (Ver. II)

----- A General Protocol for Transfecting Mammalian Cell

- 100 µl
- 500 µl
- 1000 µl



9601 Medical Center Drive
Rockville MD 20850
FAX. 301-560-4919
TEL. 301-330-5966
Toll Free. 1-(866)-918-6812
Email: info@signagen.com
Web: www.signagen.com

This product is for laboratory research ONLY and not for diagnostic use

Introduction:

Based on our innovative and proprietary lipid-conjugation technology, PowerFect™ Transfection Kit is a liposome based DNA & siRNA delivery tool which was formulated with our proprietary pH Dependent Conformational Change (PDCC) technology to give efficient and reproducible gene knockdown on variety of mammalian cells. PowerFect™ Transfection Kit is the most powerful yet very gentle gene delivery tool for a variety of applications including plasmid DNA and/or siRNA for most of mammalian cell types. Compared with leading products in the market, PowerFect™ is more cost-effective and always provides higher transfection efficiency with less cytotoxicity.

Contents Per Kit:

- 1x 1.0 mL of PowerFect™ In Vitro Transfection Reagent
- 1x 8.0 mL of PowerFect™ Transfection Buffer (5x)

Important Guidelines for Transfection:

- PowerFect™ reagent was formulated for DNA and siRNA transfection. The following standard protocol is given for DNA and siRNA transfection to mammalian cells. For a protocol of siRNA/DNA co-transfection, please email us at info@signagen.com
- For better efficiency, choosing PowerFect™ Transfection Buffer working solution (1x) is a must.
- To lower cytotoxicity, transfect cells in presence of serum (10%) and antibiotics.

Part I. A General Protocol for DNA Transfection.

Step I. Preparation of Working Solution of PowerFect™ Transfection Buffer (1x)

PowerFect™ Transfection Buffer (5x) is provided as 5x concentrated stock solution. To make working solution, dilute one part of the stock solution with 4 parts of ddH₂O. The PowerFect™ Transfection Buffer (1x) working solution is stable at RT for 24 months.

Note: Always keep PowerFect™ Transfection Buffer (5x) at RT. If refrigerated, white precipitates may appear. It won't affect the transfection efficiency. After dilution with 4 parts of ddH₂O to make PowerFect™ Transfection Buffer (1x) working solution, the white precipitates will disappear. Always keep PowerFect™ Transfection Buffer working solution (1x) at RT.

Step II. Cell Seeding:

Cells should be plated 18 to 24 hours prior to transfection so that the monolayer cell density reaches to the optimal ~70% confluency at the time of transfection. Complete culture medium with serum and antibiotics is freshly added to each well 30-60 min before transfection.

Note: High serum levels (>5%) with antibiotics do NOT have inhibitory effect on transfection efficiency. We recommend using complete serum/antibiotics-containing medium to grow the cells during transfection.

Step III. Preparation of PowerFect™ -DNA Complex and Transfection Procedures:

For different cell types, the optimal ratio of PowerFect™ (µL):DNA (µg) varies from 1:1 to 3:1. We recommend using PowerFect™ (µL):DNA (µg) at 2:1 at a starting point.

The following protocol is given for transfection in 24-well plates, refer to [Table 1](#) for transfection in other culture formats.

- For each well, dilute 0.5 µg of DNA into 50 µl of PowerFect™ Transfection Buffer working solution (1x) prepared from [Step I](#). Mix by vortexing.
- Add 1.0 µl of PowerFect™ reagent, vortex briefly to mix.
- Incubate for ~10 min at RT to allow PowerFect™/DNA complex to form.

- Note:** Never keep the PowerFect™/DNA complex longer than 20 min.
- Add the PowerFect™/DNA transfection mix to the cells in serum containing medium drop wise.
 - Swirl plate gently to homogenize.
 - Check transfection efficiency 24 to 48 hours post transfection. 48 hours usually give better efficiency.

Table 1. Recommended Amounts for Different Culture Vessel Formats

Culture Dish	Culture Medium (mL)	Plasmid DNA (µg)	PowerFect™ Transfection Buffer (1x) (µL)	PowerFect™ Reagent (µL)
96-well	0.1	0.1	5	0.2
48-well	0.25	0.25	25	0.5
24-well	0.5	0.5	50	1
6-well	2	2.0	200	4
35 mm dish	2	2.0	200	4
60 mm dish	4	4.0	400	8
10 cm / T75	10	10	800	20
15 cm / T175	20	20	1600	40

Storage: PowerFect™ siRNA Transfection Reagent is stable for up to 12 months at 4 °C. Always keep PowerFect™ Transfection Buffer (5x) at RT. This item is shipped at ambient temperature.

PowerFect™ In Vitro DNA and siRNA Transfection Kit

----- A General Protocol for Transfecting Mammalian Cell

- 100 µl
- 500 µl
- 1000 µl



9601 Medical Center Drive
 Rockville MD 20850
 FAX. 301-560-4919
 TEL. 301-330-5966
 Toll Free. 1-(866)-918-6812
 Email: info@signagen.com
 Web: www.signagen.com

This product is for laboratory research ONLY and not for diagnostic use

Part II. A General Protocol for siRNA Transfection.

Step I. Preparation of Working Solution of PowerFect™ Transfection Buffer

PowerFect™ Transfection Buffer (5x) is provided as 5x concentrated stock solution. To make working solution, dilute one part of the stock solution with 4 parts of ddH₂O. The 1x PowerFect™ Transfection Buffer is stable at RT for 24 months.

Note: Always keep PowerFect™ Transfection Buffer (5x) at RT. If refrigerated, white precipitates may appear. It won't affect the transfection efficiency. After dilution with 4 parts of ddH₂O to make PowerFect™ Transfection Buffer (1x) working solution, the white precipitates will disappear. Always keep PowerFect™ Transfection Buffer working solution (1x) at RT.

Step II. Cell Seeding:

Cells should be plated 18 to 24 hours prior to transfection so that the monolayer cell density reaches to the optimal ~50% confluency at time of transfection. Complete culture medium with serum and antibiotics is freshly added to each well 30-60 min before transfection.

Note: High serum levels (>5%) with antibiotics usually do not have inhibitory effect on transfection efficiency. We recommend using complete serum/antibiotics-containing medium to grow the cells during transfection.

Step III. Preparation of PowerFect™-siRNA Complex and Transfection Procedures:

For optimal siRNA-mediated silencing, we recommend using 10-80 nM siRNA (final concentration). The following protocol is given for transfection in 6-well plate, refer to **Table 2** for transfection in other culture formats.

- For each well, dilute 20 - 160 pmoles siRNA (for a final concentration of 10 to 80 nM per well) into 200 µl of PowerFect™ Transfection Buffer (1x) prepared from **Step I**. Mix gently.
- Add 4 µl of PowerFect™ reagent, vortex briefly to mix.
- Incubate for ~10 min at RT to allow PowerFect™/siRNA complexes to form.

Note: Never keep the PowerFect™/siRNA complex longer than 20 min.

- Add the PowerFect™/siRNA transfection mix to the cells in serum-containing medium drop wise.
- Swirl plate gently to homogenize.
- Check siRNA silencing efficiency 24 to 72 hours post transfection. 48-72 hours usually give better efficiency.

Table 2. Recommended Amounts for Different Culture Vessel Formats

Culture Dish	Culture Medium (mL)	siRNA (pmoles) 10-80 nM	PowerFect™ Transfection Buffer (1x) (µL)	PowerFect™ Reagent (µL)
96-well	0.1	1 - 8	10	0.4
48-well	0.25	2.5 - 20	25	1
24-well	0.5	5 - 40	50	2
6-well	2	20 - 160	200	4
35 mm dish	2	20 - 160	200	4
60 mm dish	4	40 - 320	400	8
10 cm / T75	10	100 - 800	800	20

Storage: PowerFect™ siRNA Transfection Reagent is stable for up to 12 months at 4 °C. Always keep PowerFect™ Transfection Buffer (5x) at RT. This item is shipped at ambient temperature.