

GenMute™ Reagent Short Protocol – siRNA Transfection

DAY 1: Prepare GenMute™ Transfection Buffer (1x) working solution and cell seeding

- Dilute one part of stock solution of GenMute™ Transfection Buffer (5x) with 4 parts of ddH₂O to make GenMute™ Transfection Buffer working solution (1x).
- Cell seeding = ~50% confluency at the day of transfection.

DAY 2: Transfection mix = GenMute™ Reagent Y (μL) with X pmol siRNA at final 5.0 nM


- Perform transfection in the presence of serum/antibiotics.
- Use GenMute™ Transfection Buffer (1x) only.
- Perform siRNA transfection referring to the picture and table below.

1. Dilute **X** pmol of siRNA in **W** μL of GenMute™ Buffer (1x). Vortex briefly to mix

2. Add **Y** μL of GenMute™ reagent. Vortex to mix and incubate 10 min at RT

3. Add transfection mix to the cells in serum-containing medium

4. Measure gene expression Kd 48 to 72 h post transfection



Culture Dish	Growth Medium (mL)	Transfection Buffer (1x) (μL) - W	siRNA (pmol) Final 5.0 nM - X	GenMute™ Reagent (μL) - Y
24-well	1	50	5.0	1.2
12-well	1.5	75	7.5	2.0
6-well	2	100	10	2.4
60mm	5	300	25	7.2
10 cm / Flask75	10	600	50	14

DAY 3~4: Measure knock-down of the gene expression

Visit the following link for complete protocol

<https://signagen.com/products/genmute-sirna-transfection-reagent>

