

GenJet™ In Vitro DNA Transfection Reagent (Ver. II)

----- A Protocol for generation of Lentivirus from 293T cell

- 100 µl
- 500 µl
- 1000 µl



9601 Medical Center drive
A/R Bldg, Suite 341
Rockville MD 20850
FAX. 301-560-4919
TEL. 301-330-5966
Toll Free. 1-(866)-918-6812
Email: info@signagen.com
Web: www.signagen.com

This product is for laboratory research ONLY and not for diagnostic use

Introduction:

GenJet™ In Vitro DNA Transfection Reagent (Ver. II) is upgraded version of GenJet™ In Vitro DNA Transfection Reagent. With a new chemistry, more DNA condensing groups were released in the new version compared with old version GenJet™, leading to 3-20 times more efficient in DAN delivery. GenJet™ (Ver. II) was shown to generate lentivirus with extremely high titers from 293T cells.

Important Transfection Guidelines for New Version:

- Do NOT follow transfection procedures for GenJet old version. Read protocol for new version carefully before transfection
- For high efficiency, transfect cells at high density. 90% cell confluency is highly recommended
- To lower cytotoxicity, transfect cells in presence of serum (10%) and antibiotics
- Change medium with serum (10%) and antibiotics 5 hours post transfection is optional

Procedures for Transfecting 293T Cells:

Cell Seeding (see Table 1):

Cells should be plated 18 to 24 hours prior to transfection so that the monolayer cell density reaches to the optimal ~90% cell confluency at the time of transfection. Complete culture medium with serum and antibiotics is freshly added to each well 30-60 minutes before transfection.

Note: High serum levels (>5%) with antibiotics usually do not have inhibitory effect on transfection efficiency. For some specific 293 cells, maximal transfection efficiencies are observed in the presence of serum and antibiotics. We recommend using complete serum/antibiotics-containing medium initially.

Table 1. A Guideline for Seeding Adherent Cells Prior to Transfection in Different Culture Formats.

Culture Dishes	Surface Area (cm ²)	Number of Cells to Seed
T75 Flask	75	3.0 - 6.0 × 10 ⁶
100 mm Dish	58	2.2 - 4.4 × 10 ⁶
60 mm Dish	21	0.9 - 1.8 × 10 ⁶
35 mm Dish	9.6	3.5 - 7.0 × 10 ⁵
6-well Plate	9.6	4.0 - 8.0 × 10 ⁵
12-well Plate	3.5	1.5 - 3.0 × 10 ⁵
24-well Plate	1.9	0.8 - 1.6 × 10 ⁵
48-well Plate	1.0	4.0 - 8.0 × 10 ⁴
96-well Plate	0.3	1.2 - 2.4 × 10 ⁴

Preparation of GenJet™-DNA Complex and Transfection Procedures

The following protocol is given for transfection in 10 cm dish. For other culture formats, scale up or down per culture dish's surface. The optimal transfection conditions are given in the standard protocol described below.

- Cell confluency should be ~90 % at the day of transfection.
- For each 10 cm dish, dilute total 15 µg of DNA (7.5 µg lentivector plasmid plus 7.5 µg lentivirus packaging mix) into 500 µL of serum-free DMEM with High Glucose. Vortex gently to mix.
- For each 10 cm dish, dilute 40 µL of GenJet™ Reagent into 500 µL of serum-free DMEM with High Glucose. Vortex gently to mix.

Note: Never use Opti-MEM to dilute DNA and GenJet™ Reagent because it will disrupt transfection complex.

- Add the diluted GenJet™ Reagent immediately to the diluted DNA solution all at once. (Important: do not mix the solutions in the reverse order !)
- Vortex to mix the solution immediately followed by incubation of 10 minutes at room temperature to allow GenJet™/DNA complexes to form.
- Note:** Never keep the DNA/GenJet™ complex longer than 20 minutes
- Add the 1000 µL GenJet™/DNA complex drop-wise onto the medium in each dish and homogenize the mixture by gently swirling the plate.
- Check transfection efficiency and harvest lentivirus supernatant 24 to 72 hours post transfection.

Storage: GenJet™ DNA In Vitro Transfection Reagent (Ver. II) is stable for up to 12 months at +4 °C. This item shipped at ambient temperature